Infusion Sets - LatexFree



CODAN US Corporation supplies high quality intravenous delivery systems and components to medical facilities throughout North America.

CODAN engineering procedures assure superior product design, safety, reliability and performance. Effective and personal customer service is a top priority at CODAN, and our entire team is committed to the needs and special requirements of each client.



LatexFree Infusion Sets

CODAN LatexFree Infusion Sets are:

- Non-DEHP
- Manufactured under ISO Standards
- Usually shipped within one business day of receipt of order

CODAN Care



Special Product Features:

- Manufactured with tubus or Y-injection site
- All sets capped for safety and security
- Options available with vented or non-vented spikes
- Drip chambers available with or without filters

CODAN US Corporation –

Your Decisive Connection to Quality IV Therapy Delivery Systems

Infusion Sets - LatexFree:

Length	Tubus	Single Y-Site	Double Y-Sites	Triple Y-Sites
75"-90" Vented	A405 C400	A410	A418 A463	
75"-112" Non-Vented	B400 B415	B411	B420	B464

Product Examples:

A405

Vented IV Set, 20 drops/ml, 15 micron fluid filter, roller clamp, tubus, luer-lock, length 80", approximate priming volume 18 ml, 50/cs

B411

Non-Vented IV Set, 20 drops/ml, roller clamp, Y-injection site 10" above luer-lock, length 82", approximate priming volume 18 ml, 50/cs

CODAN offers a wide range of IV delivery products, such as:

- Infusion Sets (NeedleFree, LatexFree and StraightLine)
- Filters and Filter Sets
- Blood Administration Sets
- Extension Sets (RegularBore, MicroBore, MiniBore and Gas Sampling)
- Light-Safe® Sets
- Transfer Sets
- IV Accessories
- Mixing and Withdrawal Products
- Custom and OEM Products

For detailed information about our full line of products and services, please call us toll-free at **(800) 33-CODAN** (26326) or visit our website at **www.codanuscorp.com**.

